# TERMINAL BLOCKS

## ASI TECHNICAL GUIDELINE FOR CHANGING TERMINAL BLOCKS

65-0235 | Rev A1 | 07/17/23 | Page 1 of 2

### PURPOSE

This technical guideline demonstrates how to remove or add wires to the terminal block.


### TOOLS REQUIRED

- Small flat head screwdriver

#### IMPORTANT!

Before commencing on the following instructions, please read and follow all applicable warnings/cautions listed at the end of this technical guideline.

### INSTRUCTIONS


The square port on the Terminal Block is the spring lock for the associated round wire port.

To open the spring lock, insert a small flat head screwdriver into the square port, and press down to open the lock; maintain downward pressure to keep lock open.

Insert or remove the associated wire.

To close, remove the screwdriver and the spring lock closes.

# TERMINAL BLOCKS

## ASI Technical Guideline for Changing Terminal Blocks

65-0235 | Rev A1 | 07/17/23 | Page 2 of 2

### WARNINGS/CAUTIONS

In addition to observing the normal precautions associated with standard dental practices and procedures, the following additional precautions should be strictly noted and observed during the set-up, operation, and maintenance of this system.

#### WARNING

##### QUALIFIED PERSONNEL ONLY

The product should only be operated by qualified personnel only. The operator bears responsibility for the correct settings and proper use of the system. ASI Dental (ASI) cannot be held liable for any malfunction of this product, or performance failure and/or its designed or desired utility, nor can ASI be held liable for injuries to persons or animals, in any case when the device is misused or not operated, applied or maintained in strict accordance with user/owner instructions set out in the operation manual. In the event of any doubt or question, the user is to contact ASI for clarification or assistance.

Improperly maintained or operated systems or instruments may void the associated warranties.

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#### WARNING

##### ELECTRICAL VOLTAGE

This system is powered by high voltage electricity. Like any other electrically powered device, if it is not used properly, it can cause electrical shock. Always plug the power cord into an electrical outlet with adequate fuse protection and proper grounding. In the event of a short circuit, grounding reduces the risk of shock by providing an escape wire for the electric current. Improper grounding of the unit can result in a risk of electric shock. Always unplug the unit before doing any service or repair to the unit.

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